02.ISBI19 notebook

January 6, 2021

1 Second exercice: iid Variable Density Sampling

This second notebook is intended to demonstrate the interest of using *iid variable density sampling* either from: * **optimal distributions** from the CS theory on orthonormal systems for Shannon wavelets derived in: - Chauffert et al, "Variable density compressed sensing in MRI. Theoretical vs heuristic sampling strategies", Proc. 10th IEEE ISBI 2013: 298-301 - Chauffert et al, "Variable Density Sampling with Continuous Trajectories" SIAM Imaging Sci, 2014;7(4):1992-1992)

• or from handcrafted densities parameterized by the decay η :


$$p(k_x, k_y) = 1/(k_x^2 + k_y^2)^{\eta/2}, \quad \eta \simeq 3.$$

• Author: Philippe Ciuciu (philippe.ciuciu@cea.fr)

• Date: 04/02/2019

- Target: ISBI'19 tutorial on Recent advances in acquisition and reconstruction for Compressed Sensing MRI
- \bullet Revision: 01/05/2020 for ATSI MSc hands-on session at Paris-Saclay University.


```
FOV = 0.2 \# field of view in m
pixelSize = FOV/img_size
#load data file corresponding to the target resolution
filename = "BrainPhantom%s.png" % img_size
mri_filename = op.join(dirimg_2d, filename)
mri_img = io.imread(mri_filename)
# mri_img = io.imread(mri_filename, as_gray=True)
# print(mri_img.dtype)
plt.figure()
plt.title("Brain Phantom, size = "+ str(img_size))
if mri_img.ndim == 2:
    plt.imshow(mri_img, cmap=plt.cm.gray)
else:
    plt.imshow(mri_img)
plt.title("Original brain image")
plt.show()
```




```
[16]: # Load target sampling distribution (precalculated in Matlab)
import numpy as np
from scipy.io import loadmat
img_size = 512
```

```
# see Chauffert et al, IEEE ISBI 2013 for the computation of optimal sampling
\rightarrow densities
densities = loadmat(op.join(dirimg_2d, "2d_sampling_densities.mat"))
if img_size == 512:
    opt density = densities['distrib2d N512 sym10'] # generated using
→orthogonal Symmlet10 wavelet transform for img size 512x512
else:
    opt_density = densities['distrib2d N256_sym10'] # generated using_
→orthogonal Symmlet10 wavelet transform for img size 256x256
# Generate Cartesian variable density mask
# change the value below if you want to change the final subsampling mask
threshold = 10. * opt_density.min() # sys.float_info.epsilon \simeq 2e-16
kspace_mask = np.zeros((img_size, img_size), dtype="float64")
kspace_mask = np.where(opt_density > threshold, 1, kspace_mask)
#plt.figure()
#plt.title("Optimal variable density for Shannon wavelets")
#plt.imshow(opt_density)
#plt.show()
fig, axs = plt.subplots(1, 2, figsize=(7, 4))
axs[0].imshow(opt_density, cmap='Greys_r')
axs[0].set title("Optimal variable density\nfor Shannon wavelets")
axs[1].imshow(kspace_mask, cmap='Greys_r')
axs[1].set title("Variable density sampling mask")
```

[16]: Text(0.5, 1.0, 'Variable density sampling mask')



```
[17]: #import numpy.fft as fft
      norm = "ortho"
      \#norm = None
      def fft(x):
          return np.fft.fft2(x, norm=norm)
      def ifft(x):
          return np.fft.ifft2(x, norm=norm)
      # Generate the kspace data: first Fourier transform the image
      kspace_data = np.fft.fftshift(fft(mri_img))
      #add Gaussian complex-valued random noise
      signoise = 10
      kspace_data += np.random.randn(*mri_img.shape) * signoise * (1+1j)
      # Mask data to perform subsampling
      kspace_data *= kspace_mask
      # Zero order solution
      image_rec0 = ifft(np.fft.ifftshift(kspace_data))
      fig, axs = plt.subplots(2, 2, figsize=(7, 7))
      axs[0, 0].imshow(mri img, cmap='Greys r')
      axs[0, 0].set_title("True image")
      axs[0, 1].imshow(kspace_mask, cmap='Greys_r')
      axs[0, 1].set_title("Sampling mask")
      axs[1,0].imshow(np.abs(kspace_data), cmap='gray', vmax=0.01*np.
      →abs(kspace_data).max())
      # axs[1].imshow(np.abs(np.fft.ifftshift(kspace_data)), cmap='Greys_r')
      axs[1, 0].set_title("k-space noisy data")
      axs[1, 1].imshow(np.abs(image_rec0), cmap='Greys_r')
      axs[1, 1].set_title("Zero-order recon")
      plt.show()
```



```
#print(p_decay.min())
      # change the value below if you want to change the final subsampling mask
      threshold = 2* opt_density.min() # sys.float_info.epsilon \simeq 2e-16
      kspace_mask = np.zeros((img_size,img_size), dtype="float64")
      kspace_mask = np.where(p_decay > threshold, 1, kspace_mask)
      fig, axs = plt.subplots(1, 2, figsize=(7, 4))
      axs[0].imshow(p_decay, cmap='Greys_r', vmax=0.001 * np.abs(p_decay).max())
      axs[0].set title("Handcrafted VD")
      axs[1].imshow(kspace mask, cmap='Greys r')
      axs[1].set_title("VD sampling mask")
     [[0.22507908 0.22463904 0.22419987 ... 0.22419987 0.22463904 0.22507908]
      [0.22463904 0.22419814 0.22375811 ... 0.22375811 0.22419814 0.22463904]
      [0.22419987 0.22375811 0.22331721 ... 0.22331721 0.22375811 0.22419987]
      [0.22419987 0.22375811 0.22331721 ... 0.22331721 0.22375811 0.22419987]
      [0.22463904 0.22419814 0.22375811 ... 0.22375811 0.22419814 0.22463904]
      [0.22507908 0.22463904 0.22419987 ... 0.22419987 0.22463904 0.22507908]]
     <IPython.core.display.Javascript object>
     <IPython.core.display.HTML object>
[13]: Text(0.5, 1.0, 'VD sampling mask')
[14]: # Generate the kspace data: first Fourier transform the image
      kspace_data = np.fft.fftshift(fft(mri_img))
      #add Gaussian complex-valued random noise
      signoise = 10
      kspace data += np.random.randn(*mri img.shape) * signoise * (1+1j)
      # Mask data to perform subsampling
      kspace_data *= kspace_mask
      # Zero order solution
      image rec0 = ifft(np.fft.ifftshift(kspace data))
      fig, axs = plt.subplots(2, 2, figsize=(7, 7))
      axs[0, 0].imshow(mri_img, cmap='Greys_r')
      axs[0, 0].set_title("True image")
      axs[0, 1].imshow(kspace_mask, cmap='Greys_r')
      axs[0, 1].set_title("Sampling mask")
      axs[1, 0].imshow(np.abs(kspace_data), cmap='gray', vmax=0.01*np.
      →abs(kspace data).max())
      #axs[1].imshow(np.abs(np.fft.ifftshift(kspace_data)), cmap='Greys_r')
```

```
axs[1, 0].set_title("k-space noisy data")
axs[1, 1].imshow(np.abs(image_rec0), cmap='Greys_r')
axs[1, 1].set_title("Zero-order recon")
plt.show()
```

<IPython.core.display.Javascript object>

<IPython.core.display.HTML object>